

Tabel 1. Data mentah pengamatan mikroskopis leukosit urine

Kode sampel	segera	Tunda 4 jam suhu 4°C	Tunda 4 jam + formaldehida 10%	Tunda 6 jam suhu 4°C	Tunda 6 jam + formaldehida 10%
1	1-3	0-2	0-2	0-2	0-2
2	19-21	16-18	16-18	16-18	15-17
3	5-7	3-5	3-5	1-3	1-3
4	4-6	3-5	3-5	2-4	1-3
5	12-14	10-12	11-13	7-9	8-10
6	0-2	0-2	0-2	0-2	0-2

Tabel 2. Nilai tengah hasil pengamatan mikroskopis leukosit urine

Kode sampel	segera	Tunda 4 jam suhu 4°C	Tunda 4 jam + formaldehida 10%	Tunda 6 jam suhu 4°C	Tunda 6 jam + formaldehida 10%
1	2	1	1	1	1
2	20	17	17	17	16
3	6	4	4	2	2
4	5	4	4	3	2
5	13	11	12	8	9
6	1	1	1	1	1